

Low affinity immunoglobulin gamma Fc region receptor III

Catalog No: #AP83223

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83223-1 50ug #AP83223-2 100ug #AP83223-3 1mg

Description

Product Name	Low affinity immunoglobulin gamma Fc region receptor III
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IgG Fc receptor III, Fc-gamma RIII, FcRIII, CD16
Accession No.	Uniprot ID:P27645
Uniprot	P27645
GeneID	100911825;108348047;116591;498276;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Fcgr3

Note: This product is for in vitro research use only