

## Free fatty acid receptor 2

Catalog No: #AP83224



Package Size: #AP83224-1 50ug #AP83224-2 100ug #AP83224-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                        |
|--------------------|--------------------------------------------------------|
| Product Name       | Free fatty acid receptor 2                             |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli                                                 |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Hu                                                     |
| Other Names        | G-protein coupled receptor 43,FFA2,GPCR43,GPR43        |
| Accession No.      | Uniprot ID:O15552                                      |
| Uniprot            | O15552                                                 |
| GeneID             | 2867;                                                  |
| Tag Info           | His                                                    |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FFAR2

Note: This product is for in vitro research use only