

Fibroblast growth factor 2

Catalog No: #AP83225



Package Size: #AP83225-1 50ug #AP83225-2 100ug #AP83225-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	FGF-2,Basic fibroblast growth factor,bFGF,Heparin-binding growth factor 2,HBGF-2,Fgf-2
Accession No.	Uniprot ID:P13109
Uniprot	P13109
GeneID	54250;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fgf2

Note: This product is for in vitro research use only