

Protein FAM19A5

Catalog No: #AP83227



Package Size: #AP83227-1 50ug #AP83227-2 100ug #AP83227-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein FAM19A5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Chemokine-like protein Tafa-5,Tafa5
Accession No.	Uniprot ID:Q91WE9
Uniprot	Q91WE9
GeneID	106014;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fam19a5

Note: This product is for in vitro research use only