

Ferritin heavy chain

Catalog No: #AP83230



Package Size: #AP83230-1 50ug #AP83230-2 100ug #AP83230-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ferritin heavy chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Ferritin H subunit
Accession No.	Uniprot ID:P08267
Uniprot	P08267
GeneID	395970;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FTH

Note: This product is for in vitro research use only