

Transcription factor GATA-4

Catalog No: #AP83236

Package Size: #AP83236-1 50ug #AP83236-2 100ug #AP83236-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor GATA-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DNA-binding protein GATA-GT2,GATA-binding factor 4
Accession No.	Uniprot ID:P46152
Uniprot	P46152
GeneID	54254;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gata4

Note: This product is for in vitro research use only