

Estrogen-related receptor gamma

Catalog No: #AP83237



Package Size: #AP83237-1 50ug #AP83237-2 100ug #AP83237-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Estrogen-related receptor gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ERR gamma-2,Estrogen receptor-related protein 3,Nuclear receptor subfamily 3 group B member 3,ERR3,ERRG2,KIAA0832,NR3B3
Accession No.	Uniprot ID:P62508
Uniprot	P62508
GeneID	2104;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ESRRG

Note: This product is for in vitro research use only