

LisH domain-containing protein FOPNL

Catalog No: #AP83244



Package Size: #AP83244-1 50ug #AP83244-2 100ug #AP83244-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LisH domain-containing protein FOPNL
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	FGFR1OP N-terminal-like protein,RCJMB04_21j4
Accession No.	Uniprot ID:Q5ZJ24
Uniprot	Q5ZJ24
GeneID	416598;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FOPNL

Note: This product is for in vitro research use only