

Squalene synthase

Catalog No: #AP83246



Package Size: #AP83246-1 50ug #AP83246-2 100ug #AP83246-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Squalene synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SQS,FPP:FPP farnesyltransferase,Farnesyl-diphosphate farnesyltransferase
Accession No.	Uniprot ID:P37268
Uniprot	P37268
GeneID	2222;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FDFT1

Note: This product is for in vitro research use only