

Fermitin family homolog 2

Catalog No: #AP83247



Package Size: #AP83247-1 50ug #AP83247-2 100ug #AP83247-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fermitin family homolog 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Kindlin-2,Pleckstrin homology domain-containing family C member 1,Plekhc1
Accession No.	Uniprot ID:Q8CIB5
Uniprot	Q8CIB5
GeneID	218952;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fermt2

Note: This product is for in vitro research use only