

Fibulin-2

Catalog No: #AP83252



Package Size: #AP83252-1 50ug #AP83252-2 100ug #AP83252-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fibulin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FIBL-2
Accession No.	Uniprot ID:P98095
Uniprot	P98095
GeneID	2199;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FBLN2

Note: This product is for in vitro research use only