

Flap endonuclease 1

Catalog No: #AP83253



Package Size: #AP83253-1 50ug #AP83253-2 100ug #AP83253-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Flap endonuclease 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FEN-1, DNase IV, Flap structure-specific endonuclease 1, Maturation factor 1, MF1, RAD2
Accession No.	Uniprot ID:P39748
Uniprot	P39748
GeneID	2237;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: FEN1

Note: This product is for in vitro research use only