

Fibulin-5

Catalog No: #AP83256



Package Size: #AP83256-1 50ug #AP83256-2 100ug #AP83256-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibulin-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	FIBL-5,Developmental arteries and neural crest EGF-like protein,Dance,Embryonic vascular EGF repeat-containing protein,EVEC,Dance
Accession No.	Uniprot ID:Q9WVH8
Uniprot	Q9WVH8
GeneID	29158;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fbln5

Note: This product is for in vitro research use only