

Immunoglobulin alpha Fc receptor

Catalog No: #AP83257



Package Size: #AP83257-1 50ug #AP83257-2 100ug #AP83257-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Immunoglobulin alpha Fc receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IgA Fc receptor,CD89,CD89
Accession No.	Uniprot ID:P24071
Uniprot	P24071
GeneID	2204;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FCAR

Note: This product is for in vitro research use only