

Histone H2AX

Catalog No: #AP83259



Package Size: #AP83259-1 50ug #AP83259-2 100ug #AP83259-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Histone H2AX
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	H2a/x,Histone H2A.X,H2AX
Accession No.	Uniprot ID:P16104
Uniprot	P16104
GeneID	3014;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:H2AFX

Note: This product is for in vitro research use only