

Oxygen-dependent coproporphyrinogen-III oxidase, mitochondrial

Catalog No: #AP83260

Package Size: #AP83260-1 50ug #AP83260-2 100ug #AP83260-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxygen-dependent coproporphyrinogen-III oxidase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	COX,Cpo
Accession No.	Uniprot ID:Q3B7D0
Uniprot	Q3B7D0
GeneID	304024;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cpox

Note: This product is for in vitro research use only