

Filamin-C

Catalog No: #AP83261



Package Size: #AP83261-1 50ug #AP83261-2 100ug #AP83261-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Filamin-C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FLN-C,ABP-280-like protein,ABP-L,Actin-binding-like protein,Filamin-2,Gamma-filamin,ABPL,FLN2
Accession No.	Uniprot ID:Q14315
Uniprot	Q14315
GeneID	2318;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FLNC

Note: This product is for in vitro research use only