

Far upstream element-binding protein 1

Catalog No: #AP83263



Package Size: #AP83263-1 50ug #AP83263-2 100ug #AP83263-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Far upstream element-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FBP,DNA helicase V,hDH V
Accession No.	Uniprot ID:Q96AE4
Uniprot	Q96AE4
GeneID	8880;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FUBP1

Note: This product is for in vitro research use only