

Endothelin-3

Catalog No: #AP83264



Package Size: #AP83264-1 50ug #AP83264-2 100ug #AP83264-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Endothelin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	ET-3,Preproendothelin-3,PPET3
Accession No.	Uniprot ID:Q765Z4
Uniprot	Q765Z4
GeneID	403406;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EDN3

Note: This product is for in vitro research use only