

Early growth response protein 1

Catalog No: #AP83267



Package Size: #AP83267-1 50ug #AP83267-2 100ug #AP83267-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Early growth response protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EGR-1,AT225,Nerve growth factor-induced protein A,NGFI-A,Transcription factor ETR103,Transcription factor Zif268,Zinc finger protein 225,Zinc finger protein Krox-24,KROX24,ZNF225
Accession No.	Uniprot ID:P18146
Uniprot	P18146
GeneID	1958;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:EGR1

Note: This product is for in vitro research use only