

## Epithelial membrane protein 1

Catalog No: #AP83271



Package Size: #AP83271-1 50ug #AP83271-2 100ug #AP83271-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Epithelial membrane protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	EMP-1, Tumor-associated membrane protein, Tmp
Accession No.	Uniprot ID:P47801
Uniprot	P47801
GeneID	13730;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Emp1

Note: This product is for in vitro research use only