

Receptor tyrosine-protein kinase erbB-2

Catalog No: #AP83273



Package Size: #AP83273-1 50ug #AP83273-2 100ug #AP83273-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor tyrosine-protein kinase erbB-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Epidermal growth factor receptor-related protein,Proto-oncogene Neu,Proto-oncogene c-ErbB-2,p185erbB2,p185neu,CD340,Neu
Accession No.	Uniprot ID:P06494
Uniprot	P06494
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ErbB2

Note: This product is for in vitro research use only