

Peptidyl-prolyl cis-trans isomerase FKBP5

Catalog No: #AP83275



Package Size: #AP83275-1 50ug #AP83275-2 100ug #AP83275-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peptidyl-prolyl cis-trans isomerase FKBP5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PPlase FKBP5,51 kDa FK506-binding protein,51 kDa FKBP,54 kDa progesterone receptor-associated immunophilin,Androgen-regulated protein 6,FF1 antigen,FK506-binding protein 5,FKBP-5,FKBP54,p54,HSP90-binding immunophilin,Rotamase,AIG6,FKBP51
Accession No.	Uniprot ID:Q13451
Uniprot	Q13451
GeneID	2289;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FKBP5

Note: This product is for in vitro research use only