

Forkhead box protein O3

Catalog No: #AP83276



Package Size: #AP83276-1 50ug #AP83276-2 100ug #AP83276-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Forkhead box protein O3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AF6q21 protein,Forkhead in rhabdomyosarcoma-like 1,FKHRL1,FOXO3A
Accession No.	Uniprot ID:O43524
Uniprot	O43524
GeneID	2309;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FOXO3

Note: This product is for in vitro research use only