

# Ficolin-2

Catalog No: #AP83280



Package Size: #AP83280-1 50ug #AP83280-2 100ug #AP83280-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ficolin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Collagen/fibrinogen domain-containing protein 2,Ficolin-B,Ficolin-beta,L-ficolin
Accession No.	Uniprot ID:Q5I2E5
Uniprot	Q5I2E5
GeneID	497016;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FCN2

Note: This product is for in vitro research use only