

Fibroleukin

Catalog No: #AP83282



Package Size: #AP83282-1 50ug #AP83282-2 100ug #AP83282-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroleukin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cytotoxic T-lymphocyte-specific protein,Fibrinogen-like protein 2,Prothrombinase,Fiblp
Accession No.	Uniprot ID:P12804
Uniprot	P12804
GeneID	14190;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fgl2

Note: This product is for in vitro research use only