

Gastric inhibitory polypeptide receptor

Catalog No: #AP83285



Package Size: #AP83285-1 50ug #AP83285-2 100ug #AP83285-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gastric inhibitory polypeptide receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GIP-R, Glucose-dependent insulinotropic polypeptide receptor
Accession No.	Uniprot ID:P43219
Uniprot	P43219
GeneID	25024;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Gipr

Note: This product is for in vitro research use only