

Free fatty acid receptor 4

Catalog No: #AP83286



Package Size: #AP83286-1 50ug #AP83286-2 100ug #AP83286-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Free fatty acid receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	G-protein coupled receptor 120, Omega-3 fatty acid receptor 1, Gpr120, O3far1
Accession No.	Uniprot ID:Q2AC31
Uniprot	Q2AC31
GeneID	294075;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ffar4

Note: This product is for in vitro research use only