

Peptidyl-prolyl cis-trans isomerase FKBP5

Catalog No: #AP83288



Package Size: #AP83288-1 50ug #AP83288-2 100ug #AP83288-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peptidyl-prolyl cis-trans isomerase FKBP5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PPIase FKBP5,51 kDa FK506-binding protein,51 kDa FKBP,FK506-binding protein 5,FKBP-5,Rotamase,Fkbp51
Accession No.	Uniprot ID:Q64378
Uniprot	Q64378
GeneID	14229;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fkbp5

Note: This product is for in vitro research use only