

Transcription factor HES-1

Catalog No: #AP83293



Package Size: #AP83293-1 50ug #AP83293-2 100ug #AP83293-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor HES-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Hairy and enhancer of split 1
Accession No.	Uniprot ID:Q3ZBG4
Uniprot	Q3ZBG4
GeneID	539547;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HES1

Note: This product is for in vitro research use only