

# H-2 class II histocompatibility antigen gamma chain

Catalog No: #AP83294



Package Size: #AP83294-1 50ug #AP83294-2 100ug #AP83294-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	H-2 class II histocompatibility antigen gamma chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ia antigen-associated invariant chain, Ii, MHC class II-associated invariant chain, CD74, Ii
Accession No.	Uniprot ID:P04441
Uniprot	P04441
GeneID	16149;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Cd74

Note: This product is for in vitro research use only