

Receptor tyrosine-protein kinase erbB-2

Catalog No: #AP83296



Package Size: #AP83296-1 50ug #AP83296-2 100ug #AP83296-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Receptor tyrosine-protein kinase erbB-2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Metastatic lymph node gene 19 protein,MLN 19,Proto-oncogene Neu,Proto-oncogene c-ErbB-2,Tyrosine kinase-type cell surface receptor HER2,p185erbB2,CD340,HER2,MLN19,NEU,NGL |
| Accession No. | Uniprot ID:P04626 |
| Uniprot | P04626 |
| GeneID | 2064; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ERBB2

Note: This product is for in vitro research use only