

Carboxylesterase 5A

Catalog No: #AP83297



Package Size: #AP83297-1 50ug #AP83297-2 100ug #AP83297-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carboxylesterase 5A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Carboxylesterase-like urinary excreted protein homolog,Cauxin,CES7
Accession No.	Uniprot ID:Q6AW47
Uniprot	Q6AW47
GeneID	445456;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CES5A

Note: This product is for in vitro research use only