

Alpha-2-HS-glycoprotein

Catalog No: #AP83298



Package Size: #AP83298-1 50ug #AP83298-2 100ug #AP83298-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpha-2-HS-glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Fetuin-A,FETUA
Accession No.	Uniprot ID:P29700
Uniprot	P29700
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:AHSG

Note: This product is for in vitro research use only