

Free fatty acid receptor 4

Catalog No: #AP83300



Package Size: #AP83300-1 50ug #AP83300-2 100ug #AP83300-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Free fatty acid receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	G-protein coupled receptor 120,G-protein coupled receptor GT01,Omega-3 fatty acid receptor 1,Gpr120,O3far1
Accession No.	Uniprot ID:Q7TMA4
Uniprot	Q7TMA4
GeneID	107221;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ffar4

Note: This product is for in vitro research use only