

Fibroblast growth factor 6

Catalog No: #AP83301



Package Size: #AP83301-1 50ug #AP83301-2 100ug #AP83301-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Fibroblast growth factor 6 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | FGF-6,Heparin secretory-transforming protein 2,HST-2,Heparin-binding growth factor 6,HBGF-6,HST2,HSTF2 |
| Accession No. | Uniprot ID:P10767 |
| Uniprot | P10767 |
| GeneID | 2251; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGF6

Note: This product is for in vitro research use only