

# Glutathione S-transferase P

Catalog No: #AP83302



Package Size: #AP83302-1 50ug #AP83302-2 100ug #AP83302-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glutathione S-transferase P
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	GST P1-1,GST class-pi
Accession No.	Uniprot ID:P80031
Uniprot	P80031
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:GSTP1

Note: This product is for in vitro research use only