

Fibroblast growth factor-binding protein 1

Catalog No: #AP83304



Package Size: #AP83304-1 50ug #AP83304-2 100ug #AP83304-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	FGF-BP
Accession No.	Uniprot ID:O70514
Uniprot	O70514
GeneID	14181;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Fgfbp1

Note: This product is for in vitro research use only