

Metabotropic glutamate receptor 5

Catalog No: #AP83305



Package Size: #AP83305-1 50ug #AP83305-2 100ug #AP83305-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Metabotropic glutamate receptor 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	mGluR5,Gprc1e,Mglur5
Accession No.	Uniprot ID:Q3UVX5
Uniprot	Q3UVX5
GeneID	108071;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Grm5

Note: This product is for in vitro research use only