

Solute carrier family 2, facilitated glucose transporter member 1

Catalog No: #AP83307

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83307-1 50ug #AP83307-2 100ug #AP83307-3 1mg

Description

Product Name	Solute carrier family 2, facilitated glucose transporter member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Glucose transporter type 1, erythrocyte/brain, GLUT-1, GLUT-I, GLUT1
Accession No.	Uniprot ID:P27674
Uniprot	P27674
GeneID	282356;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC2A1

Note: This product is for in vitro research use only