

Glycerol kinase

Catalog No: #AP83312



Package Size: #AP83312-1 50ug #AP83312-2 100ug #AP83312-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glycerol kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	GK,ATP:glycerol 3-phosphotransferase
Accession No.	Uniprot ID:Q0IID9
Uniprot	Q0IID9
GeneID	505987;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GK

Note: This product is for in vitro research use only