

Gastrotropin

Catalog No: #AP83313



Package Size: #AP83313-1 50ug #AP83313-2 100ug #AP83313-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gastrotropin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	GT,Fatty acid-binding protein 6,ileal lipid-binding protein,ILBP,Porcine ileal peptide,PIP,ILBP,ILLBP
Accession No.	Uniprot ID:P10289
Uniprot	P10289
GeneID	397423;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FABP6

Note: This product is for in vitro research use only