

Frizzled-9

Catalog No: #AP83315



Package Size: #AP83315-1 50ug #AP83315-2 100ug #AP83315-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Frizzled-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fz-9,CD349,Fzd3
Accession No.	Uniprot ID:Q9R216
Uniprot	Q9R216
GeneID	14371;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Fzd9

Note: This product is for in vitro research use only