

Trans-acting T-cell-specific transcription factor GATA-3

Catalog No: #AP83320

Package Size: #AP83320-1 50ug #AP83320-2 100ug #AP83320-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Trans-acting T-cell-specific transcription factor GATA-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	GATA-binding factor 3
Accession No.	Uniprot ID:Q08DV0
Uniprot	Q08DV0
GeneID	505169;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GATA3

Note: This product is for in vitro research use only