

## LisH domain-containing protein FOPNL

Catalog No: #AP83323



Package Size: #AP83323-1 50ug #AP83323-2 100ug #AP83323-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LisH domain-containing protein FOPNL
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	FGFR1OP N-terminal-like protein,FOP-related protein of 20 kDa,C16orf63,FOR20,PHSECRG2
Accession No.	Uniprot ID:Q96NB1
Uniprot	Q96NB1
GeneID	123811;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FOPNL

Note: This product is for in vitro research use only