

# Forkhead box protein O3

Catalog No: #AP83324



Package Size: #AP83324-1 50ug #AP83324-2 100ug #AP83324-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Forkhead box protein O3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Fkhr2,Foxo3a
Accession No.	Uniprot ID:Q9WVH4
Uniprot	Q9WVH4
GeneID	56484;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Foxo3

Note: This product is for in vitro research use only