

# Protein FAM3A

Catalog No: #AP83325



Package Size: #AP83325-1 50ug #AP83325-2 100ug #AP83325-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein FAM3A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytokine-like protein 2-19,2-19,2.19
Accession No.	Uniprot ID:P98173
Uniprot	P98173
GeneID	60343;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FAM3A

Note: This product is for in vitro research use only