

RNA-binding protein FUS

Catalog No: #AP83327



Package Size: #AP83327-1 50ug #AP83327-2 100ug #AP83327-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RNA-binding protein FUS
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	75 kDa DNA-pairing protein,Oncogene FUS,Oncogene TLS,POMp75,Translocated in liposarcoma protein,TLS
Accession No.	Uniprot ID:P35637
Uniprot	P35637
GeneID	2521;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FUS

Note: This product is for in vitro research use only