

Fibroblast growth factor-binding protein 2

Catalog No: #AP83328



Package Size: #AP83328-1 50ug #AP83328-2 100ug #AP83328-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Fibroblast growth factor-binding protein 2 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Chicken |
| Other Names | FGF-BP2 |
| Accession No. | Uniprot ID:Q802A9 |
| Uniprot | Q802A9 |
| GeneID | 395094; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGFBP2

Note: This product is for in vitro research use only