

## GDNF family receptor alpha-1

Catalog No: #AP83333



Package Size: #AP83333-1 50ug #AP83333-2 100ug #AP83333-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GDNF family receptor alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	GDNF receptor alpha-1,RET ligand 1,TGF-beta-related neurotrophic factor receptor 1,Gdnfra,Retl1,Tnrr1
Accession No.	Uniprot ID:Q62997
Uniprot	Q62997
GeneID	25454;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Gfra1

Note: This product is for in vitro research use only